

DrySyn SnowStorm Application Note

Comparing performance with a conventional ice bath.
This test was carried out independently by an Asynt customer.



Objective:

To observe the active cooling of a DrySyn SnowStorm unit compared to a passive ice bath, ascertaining which was better at controlling an exotherm.

DrySyn SnowStorm Set Up:

- 100 mL round bottomed flask, Methanol 50 mL
- DrySyn SnowStorm base and 100ml insert
- Huber Ministat 240 with temperature range of -45...200 °C
- PT100 sensor in DrySyn SnowStorm block (Monitoring Drysyn SnowStorm temperature)
- PTFE coated PT100 sensor in reaction flask for feedback control to Huber Ministat
- Magnetic stirrer, 700 rpm

Ice Bath Set Up:

- 100 mL round bottomed flask, Methanol 50 mL
- Plastic ice bath, 900 g ice and 400 mL water
- Magnetic Stirrer, 700 rpm
- Metal sheathed PT100 probe monitoring ice temperature
- Metal sheathed PT100 probe monitoring solution temperature

No insulation used.


Experiment:

Both systems were set up and monitored using a datalogger.

Both systems were equilibrated until the solution temperature was stable at 0 °C.

At time = 7.9 minutes Methanol (rt, 20 mL) was added to each reaction to simulate an exotherm. The reactions were monitored and the results recorded. (See Figure 1. overleaf)

Figure 1:



Conclusions:

- DrySyn SnowStorm controlled the high exothermic temperatures more effectively than the ice bath. The temperature in the DrySyn SnowStorm reached 13 °C whereas in the ice bath it reached 20.7 °C
- DrySyn SnowStorm rapidly cooled the temperature to 0 °C in only 5 minutes. The ice bath failed to reach 0 °C. However, after 21 minutes, the ice bath had reached 1 °C
- Figure 1 shows that when controlling an exotherm, the DrySyn SnowStorm will be more effective due to the

Further Information:

DrySyn SnowStorm ONE



Huber Ministat 240

