Calorimetry

Boiling Point



Boiling point of engine coolants Covers the determination of the equilibrium boiling point of engine coolants

Art. LT/BP-232000/M Boiling Point ASTM D1120

- 100 ml round-bottom short-neck, heat resistant glass flask with standard taper, female ground-glass joint and side-entering tube for the introduction of the thermometer
- Condenser: water-cooled reflux glass-tube type
- Heating mantle with rod, clamp and pliers
- Thermometer ASTM 2C

Spare Parts

LAB-102-321: flask LAB-102-322: condenser T-AS2C: thermometer ASTM 2C

CE

