

# Haier -40 °C Deep Freezers



These products are applicable in research institutes laboratories blood banks, hospitals epidemic prevention services and so on.

## ◆Temperature Control

- Microprocessor control,digital temperature display,adjustable temperature range: -20 ~ -40°C.
- Over/under temperature controller with adjustable alarm set points.

## ◆Safety Control

- Malfunction alarms:over/under temperature,sensor error.
- Alarm representations:buzzing,flashing.

## ◆Refrigeration System

- Danfoss compressor and Germacy EBM condensation fan to get a high cooling effect.
- Performance(can reduce to -40°C within 2 hours under the ambient temperature of 35°C).
- High-density insulating layer.
- Reasonable design of evaporator and condensor system.
- CFC free refrigerant.

## ◆Ergonomic Design

- Safety lock to prevent unauthorized access.
- Ç 40 access port for operation under 198~ 252 V.
- Broad voltage band for operation under 198~ 252V.
- Stepped rust-resisting interior cabinet design to keep various baskets

## Parameter

Model No.	DW-40W380
Voltage(V/Hz)	220/50
Input power(w)	550
Temp range(°C)	-20~ -40
Exterior Dimension(mm)(WxDxH)	1550X633X838
Interior Dimensions (mm)(WxDxH)	1376X457X625
Litres(L)	380
Net/Gross(kg)	