

SampleGenie 4 and ROCKET SYNERGY

- 250ml of sample into the vial of your choice

SampleGenie 4 is the latest version of Genevac's unique SampleGenie for vials, designed for use in the **Rocket Synergy** and **Rocket 4D Synergy** Evaporators.

SampleGenie technology enables large sample volumes to be dried directly into vials, eliminating a number of time-consuming sample handling steps and the attendant risk of errors. The new 4th generation further enhances the exceptional capabilities of the original design.

New for SampleGenie 4:

- Redesigned flask collar and new lower rotor speed methods reduce stress on the vial – the likelihood of vial breakage owing to variable vial quality is significantly reduced
- Precision tightening method, combined with new seal materials, enables you to achieve optimal tightness with complete confidence

If you haven't tried the Rocket with SampleGenie for vials before, now is the time to experience the difference in evaporation performance and operational savings that they make possible.

Go to www.genevac.com for more information.



Precision tightening method, with a visual guide, enables the user to achieve optimal tightness without risking vial breakage.



Flask with sample volume of 250ml.

Choice of one piece or two piece, Viton or Perfluoroelastomer seals for leak free operation.

Visual guide for precise tightness (see image on left).

New vial adapter insert to suit a range of vials from 12mm to 28mm dia. and up to 70mm tall.