


Zr Low Form Crucibles



Code	Capacity	A	B	C
ZLF05	5ml	22mm	22mm	14mm
ZLF10	10ml	27mm	24mm	22mm
ZLF20	20ml	36mm	33mm	28mm
ZLF25	25ml	40mm	28mm	32mm
ZLF30	30ml	38mm	35mm	32mm
ZLF35	35ml	39mm	38mm	32mm
ZLF40	40ml	41mm	39mm	36mm
ZLF45	45ml	48mm	39mm	36mm
ZLF50	50ml	50mm	39mm	38mm
ZLF55	55ml	51mm	43mm	38mm
ZLF70	70ml	57mm	50mm	53mm
ZLF100	100ml	57mm	50mm	53mm